

BEAM SPLITTERS

OP-TV

DOUBLE THE IMAGE IN THE EASIEST WAY

- Resistant
- Bright
- Focusable



OPTOMIC offers a complete range of video solutions for showing live images on a monitor, taking pictures or capturing videos, as well as the possibility of transferring the images to a PC or laptop.

The beam splitters are the perfect complement to our microscopes and colposcopes and are very easy to operate.

BEAM SPLITTERS **OP-TV**



■ OP-TV CLICK

This beam splitter has a quick connection for easy attachment to any CCD camera for endoscopy with endocoupler.

A versatile solution for combining endoscopy and microscopy/colposcopy.



■ OP-TV DC

Beam splitter with digital photographic camera. Camera with LCD screen with adjustable angle, full HD video recording and HDMI output.

The best choice for taking still pictures with the advantage of being wireless.



■ OP-TV8

Compact system with built-in camera.

Video outputs: 1 x composite video and 1 x S-video (Y/C)

An elegant option at an affordable price.



■ OP-TV10 HD

Compact system with integrated HD camera.

Video outputs: 1 x HDMI (for the monitor) and 1 x USB 3.0 (for the computer)

The professional solution with HD quality.

OPTOMIC ESPAÑA S.A.

Madroño, 4. E-28770 Colmenar Viejo, Madrid, Spain

Tel: +34 91 804 44 44 • Fax: +34 91 803 84 14

export@optomic.com • www.optomic.com



producing quality

Specifications are subject to change without notice

