



- 40x optical zoom
- Great luminosity
- Live View screen integrated into the head
- Image and video capture



ColpoScreen expands the range of colposcopy equipment manufactured by OPTOMIC by adding high-resolution digital video technology. The ColpoScreen is a compact, easy-to-use Full HD device compatible with most image and video capture, display and management systems in the market.

ColpoScreen

VIDEOCOLPOSCOPE WITH INTEGRATED SCREEN



COMPACT HEAD



LIVE VIEW SCREEN
AND USER-FRIENDLY CONTROL



INTEGRATED LED LIGHTING
SYSTEM IN THE HEAD

GENERAL CHARACTERISTICS

- Camera:
 - Full HD resolution (1080P)
 - 40x optical zoom
 - Automatic control of gain and light intensity
 - Automatic focus with manual control option
- 3.5" Live View screen integrated into the head
- High intensity LED lighting
- Light intensity control
- Green filter (LED)
- White balance control
- HDMI 1080p/30 output for monitor
- 5-wheeled rolling floor stand
- Power: 100-240 VAC/50-60Hz

OPTIONAL ACCESSORIES

- Monitor support
- 22" Monitor
- "All in one" monitor with image and video capture.

OPTOMIC ESPAÑA S.A.

Madroño 4. E-28770, Colmenar Viejo. Madrid - Spain
Tel: +34 91 804 44 44
export@optomic.com
www.optomic.com



Specifications are subject to change without prior notice

