

ENT RIGID ENDOSCOPES OTO-LARYNGOSCOPES

OUTSTANDING QUALITY IN ENDOSCOPIC VISION


- High resolution
- Clear images
- High resistance



OPTOMIC offers a complete range of rigid otoscopes and laryngoscopes for the inspection of the corresponding area. They allow for the diagnosis and serve as a method for minimal invasive surgical intervention. All the ENT OPTOMIC endoscopes are of proven quality and offer an exceptional image brightness to guarantee the satisfaction of the specialist.

ENT RIGID ENDOSCOPES


OTO-LARYNGOSCOPES



ROD LENSES SYSTEM



STAINLESS STEEL BODY



ADAPTERS



RESILIENT SAPPHIRE LENSES




PEEK EYEPIECE

| | | | 0° | 30° | 70° | 70° | 90° |
|--------------|--------|--------|----|-----|-----|-----|-----|
| | | | | | | | |
| OTOSCOPE | 1,9 mm | 100 mm | ● | ● | | | |
| | 2,7 mm | 75 mm | ● | ● | ● | | |
| | 4 mm | 50 mm | ● | ● | | | |
| LARYNGOSCOPE | 4 mm | 185 mm | | | | ● | |
| | 6 mm | 158 mm | | | | ● | |
| | 10 mm | 158 mm | | | | ● | ● |

● Non-Autoclavable ● Autoclavable

- Rod lenses system and stainless steel sealed by laser welding.
- Adapters for connection to type *Storz/Olympus*, *Wolf* and *ACMI* fiber optic cables.
- Autoclavable endoscopes with lenses made of sapphire of high resistance to scratches.



Specifications are subject to change without notice

