

TRAINING MICROSCOPES

OP-TRAINING 3 / OP-TRAINING 5

BEST CHOICE FOR PRACTICAL MICROSCOPY TRAINING

- Specially designed for training and laboratory
- Great mobility in a reduced space
- Compact and versatile



OPTOMIC[®]
producing quality

The OP-TRAINING 3 and OP-TRAINING 5 microscopes series have been designed to meet the specific requirements of laboratory and microscopy training centers. The microscopes offer excellent optical and mechanical features with a long lasting, high-power led illumination.

TRAINING MICROSCOPE **OP-TRAINING 3 / OP-TRAINING 5**



On/off and light adjustment



High power LED illumination



Inclinable binocular (optional)



Stainless steel trays (optional)

HEAD	Stereoscopic head with micrometric focus
BINOCULAR	Choose between straight, 45° inclined or 0°-240° inclinable
ILLUMINATION SYSTEM	LED integrated in the head
EYEPIECES	10x or 12,5x
MAGNIFICATION CHANGER	3-steps (OP-TRAINING 3) and 5-steps (OP-TRAINING 5)
OBJECTIVES	250mm standard (200mm optional)
FILTER	Amber or green, depending on the specialty
ARMS	Exclusive Soft-Move system providing smooth and precise mobility thanks to its high resistance conical cylinder bearings
TABLE MOUNTING	Fixation clamp (with opening of 6,5 cm) or 3 x M8 screws
POWER SUPPLY CONNECTION	100/240V 50/60Hz

OPTIONAL ACCESSORIES	Monitor Monitor support Instrument tray with rotating support Temporal bone holder Beam splitters (various options) Digitally Soft image and video capture software
----------------------	--

The **OP-TRAINING 3** and **OP-TRAINING 5** microscopes are essential in acquiring skills and gaining competence through practice.



OPTOMIC ESPAÑA S.A.

Madroño, 4. E-28770 Colmenar Viejo, Madrid, Spain
Tel: +34 91 804 44 44 • Fax: +34 91 803 84 14
export@optomic.com • www.optomic.com

