

# LED HEADLIGHT

## OP-HL15

ACCURATE LIGHT INTENSITY AND  
COLOR REPRODUCTION FOR EVERY OCCASION

- Wireless
- Comfortable
- High intensity  
light



**OPTOMIC**  
producing quality

Lightweight and unobtrusive headlight that uses powerful LED technology to produce intense white light, while maintaining color fidelity. The headlight does not require cables, making it very comfortable to use.

## LED HEADLIGHT **OP-HL15**



LIGHT AND SLEEK DESIGN



BATTERY HOLDER



LIGHT INTENSITY ADJUSTMENT

- Operating time of 9 hours.
- Pure white light of high intensity.
- Light intensity at output: 850,000 lux.
- Light intensity at 150mm: 100,000 lux.
- Vertical adjustment of illumination angle.
- Rechargeable lithium-ion battery.
- Weight: 258 g (including battery).
- Adjustable leather headband.

### **OPTOMIC ESPAÑA S.A.**

Madroño, 4. E-28770 Colmenar Viejo. Madrid, Spain  
Tel: +34 91 804 44 44  
export@optomic.com  
www.optomic.com



Specifications are subject to change without notice