

RIGID ENDOSCOPES

HYSTEROSCOPES

OUTSTANDING QUALITY IN ENDOSCOPIC VISION


- High resolution
- Clear images
- High resistance



OPTOMIC offers a complete range of rigid hysteroscopes for the inspection of the uterine cavity. They allow for the diagnosis of intrauterine pathologies and serve as a method for surgical intervention. All the OPTOMIC hysteroscopes are of proven quality and offer an exceptional image brightness to guarantee the satisfaction of the specialist.

RIGID ENDOSCOPES


HYSTEROSCOPES




ROD LENSES SYSTEM



STAINLESS STEEL BODY



ADAPTERS



RESILIENT SAPPHIRE LENSES



PEEK EYEPIECE

HYSTEROSCOPE	Ø	Length	0°	12°	30°	70°
	2,7mm	270mm	● ●		● ●	
	4mm	300mm	● ●	● ●	● ●	● ●

● Non-Autoclavables ● Autoclavables

- Endoscope with rod lenses system.
- Adapters for connection to type *Storz/Olympus*, *Wolf* and *ACMI* fiber optic cables.
- Body made of stainless steel sealed by laser welding.
- Autoclavable endoscopes with lenses made of sapphire of high resistance to scratches.
- Eyepiece of PEEK resistant to heat.

OPTOMIC ESPAÑA S.A.

Madroño, 4. E-28770 Colmenar Viejo, Madrid, Spain
Tel: +34 91 804 44 44 • Fax: +34 91 803 84 14
export@optomic.com • www.optomic.com



Specifications are subject to change without notice

