

RIGID ENDOSCOPES LAPAROSCOPES


OUTSTANDING QUALITY IN ENDOSCOPIC VISION

- High resolution
- Clear images
- High resistance



OPTOMIC offers a complete range of rigid laparoscopes for diagnostic purposes or minimal invasive surgical procedure within the abdominal or pelvic cavities. All the OPTOMIC laparoscopes are of proven quality and offer an exceptional image brightness to guarantee the satisfaction of the laparoscopic specialist.


RIGID ENDOSCOPES LAPAROSCOPES




ROD LENSES SYSTEM



STAINLESS STEEL BODY



ADAPTERS



RESILIENT SAPPHIRE LENSES



PEEK EYEPiece

	Ø	Length	0°	30°	45°
LAPAROSCOPE	5mm	300mm	● ●	● ●	● ●
	10mm	330mm	● ●	● ●	● ●

● Non-Autoclavables ● Autoclavables

- Endoscope with rod lenses system.
- Adapters for connection to type *Storz/Olympus*, *Wolf* and *ACMI* fiber optic cables.
- Body made of stainless steel sealed by laser welding.
- Autoclavable endoscopes with lenses made of sapphire of high resistance to scratches.
- Eyepiece of PEEK resistant to heat.

OPTOMIC ESPAÑA S.A.

Madroño, 4. E-28770 Colmenar Viejo, Madrid, Spain
 Tel: +34 902 22 11 77 / +34 91 803 57 14 • Fax: +34 902 22 11 78
 export@optomic.com • www.optomic.com



Specifications are subject to change without notice

