

# CYSTOSCOPE

## OP-14

### GREAT VERSATILITY

- High resolution
- Easy handling
- Working channel and irrigation



OPTOMIC OP-14 flexible endoscope is designed to facilitate passage through the urethra into the bladder to examine the mucosa for diagnostic and therapeutic procedures.

CYSTOSCOPE **OP-14**



FIBROLUX LED DUO  
LIGHT SOURCE



FIBROLUX LED HP  
LIGHT SOURCE



FIBER OPTIC CABLE



INSTRUMENTS



LEAKAGE TEST DEVICE



WORKING CHANNEL

OP-14

<b>Optical system</b>	
Field of view	110° forward
Depth of field	4-50mm
Direction of view	0°
<b>Insertion tube</b>	
Outer diameter	15.5Fr
Working length	370mm
<b>Range of tip bending</b>	
Up	210°
Down	140°
<b>Biopsy channel</b>	
Working channel inner diameter	7,0Fr

**OPTOMIC ESPAÑA S.A.**

Madroño, 4. E-28770 Colmenar Viejo, Madrid, Spain  
Tel: +34 91 804 44 44 • Fax: +34 91 803 84 14  
export@optomic.com • www.optomic.com



Specifications are subject to change without notice

