

COLPOSCOPE

OP-C2

THE BEST VISION IN COLPOSCOPY

- 3 different magnifications
- Easy to handle
- Excellent illumination



The OPTOMIC OP-C2 colposcope is the perfect equipment for diagnostic observations as well as minor routine surgeries in the doctor's practice, helping the specialist in his everyday work.

COLPOSCOPE **OP-C2**



INCLINED BINOCULAR HEAD



REMOVABLE COLD LIGHT SOURCE



OPTIONAL VIDEO SYSTEMS

- 45° inclined stereoscopic binocular head.
- Easy one-hand positioning handle.
- The arm is an auto-compensated parallelogram which can be placed in any position.
- Very soft and smooth operation of all functions.
- LED cold light source with extensive lifespan of over 50,000 hours.
- Green filter for enhance diagnostic capabilities.
- Removable cold light source can be used for other applications such as endoscopy, illumination for speculum or head lamp.
- High intensity illumination with two color temperatures (6.500 °K & 4.600 °K)
- Wide-field eyepieces with dioptric adjustment, inter pupillar adjustment and protective rubber shells.
- 3-step Galilei magnification changer and micrometric focusing system.
- Working distance of 250mm with standard objective (3 other available).
- Range of magnifications: 3.5x up to 17.5x, depending on the chosen combination of objective, binocular and eyepieces.
- Different floor and wall mounts available.
- Optional monitor supports and trays for laptop or endoscopic camera.
- Optional systems for video-colposcopy or digital photography.

OPTOMIC ESPAÑA S.A.

Madroño, 4. E-28770 Colmenar Viejo, Madrid, Spain
Tel: +34 91 804 44 44 • Fax: +34 91 803 84 14
export@optomic.com • www.optomic.com



Specifications are subject to change without notice

