

LED HEADLIGHT

OP-HL22 D

ACCURATE LIGHT INTENSITY AND
COLOR REPRODUCTION FOR EVERY OCCASION

- Custom adjustment
- Wireless
- High intensity light



OPTOMIC[®]
producing quality

Headlight with helmet regulation system, adaptable to any professional. It uses powerful LED technology to produce intense white light while maintaining color fidelity. The headlight does not require cables, making it very comfortable to use.

LED HEADLIGHT **OP-HL22 D**



LIGHT DIAPHRAGM



DUAL CRANIAL ADJUSTMENT



LIGHT INTENSITY ADJUSTMENT

- Ergonomic helmet adjustment system: two controls for perimeter and head height adjustment.
- Autonomy of 9 hours.
- High intensity pure white light.
- Light intensity at output: 850,000 lux.
- Light intensity at 150mm: 100,000 lux.
- Light spot range from 15mm to 95mm (at working distances of 400mm).
- Gradual light adjustment system.
- Vertical adjustment of illumination angle.
- Weight: 396 g (including battery).
- Rechargeable lithium-ion battery.

OPTOMIC ESPAÑA S.A.

Madroño 4. E-28770, Colmenar Viejo. Madrid - Spain
Tel: +34 91 804 44 44
export@optomic.com
www.optomic.com



producing quality

Specifications are subject to change without notice

