

DOCTOR'S STOOL

OP-B21

FUNCTIONAL STOOL FOR THE DOCTOR
WITH BACKREST AND HEIGHT ADJUSTMENT

- Ergonomic
- Resistant
- Maximum stability



producing quality

The OP-B21 is the ideal stool for any medical office thanks to its ease of movement and height adjustability providing greater comfort for the doctor.

STOOL **OP-B21**



HEIGHT LEVER



ERGONOMIC BACKREST



SWIVEL CASTORS



- Height adjustment from 550mm to 710mm
- Chrome steel structure
- Lever controlled gas spring height adjustment
- Ergonomic backrest
- Round seat of Ø 390mm
- Stable PVC base with 5 legs and easy-rolling swivel castors
- Stool weight: 10 Kg
- Available in different colours



producing quality

Specifications are subject to change without notice

