

ENDOSCOPIC FURNITURE

COMPACT TROLLEY

YOUR ENDOSCOPIC EQUIPMENT IN THE MINIMUM SPACE



- Sleek and elegant
- Light design
- Easy to handle



Compact and versatile trolley with monitor support, specially designed for your endoscopic equipment using minimum space.
The perfect complement for everyday endoscopic procedures.


ENDOSCOPIC FURNITURE **COMPACT TROLLEY**



FIBERSCOPE HOLDER



CAMERA HEAD SUPPORT



MULTIDIRECTIONAL MOVEMENT



- Adjustable tray; additional tray available.
- Multidirectional flat monitor support.
- Mounting system VESA 50x50, VESA 75x75 and VESA 100x100.
- Four wheels, two with brakes.
- Total height: 1,37m.
- Back handle for easy maneuvering.
- Optional support for flexible endoscope.

OPTOMIC ESPAÑA S.A.

Madroño, 4. E-28770 Colmenar Viejo, Madrid, Spain
Tel: +34 91 804 44 44 • Fax: +34 91 803 84 14
export@optomic.com • www.optomic.com



Specifications are subject to change without notice