

HEAD LAMP

OP-LF5

GREAT EFFICIENCY IN LIGHT TRANSMISSION

- Optimally lit area
- Easy adjustment
- Comfortable



producing quality

The OP-LF5 fiber optic head lamp is the most effective way to illuminate the working area. Its low weight and easy headband adjustment make it highly comfortable.

HEAD LAMP **OP-LF5**



FOCUSABLE OPTICS




FIBER OPTIC CABLE



CARRYING CASE

- Focusable optics.
- Ø 3,5mm x 2,3m fiber optic cable.
- Two comfortable and easily adjustable padded headbands.
- Adjustable illumination axis.
- High level brightness.
- Storz type adapter to light source (other types optional).



producing quality

Specifications are subject to change without notice

