

ENDOSCOPIC FURNITURE PROFESSIONAL TROLLEY

YOUR ENDOSCOPIC EQUIPMENT IN THE MINIMUM SPACE

- Adjustable
- Light design
- Easy to handle



OPTOMIC
producing quality

Compact and versatile trolley with monitor support, specially designed for your endoscopic equipment using minimum space.
The perfect complement for everyday endoscopic procedures.

ENDOSCOPIC FURNITURE

PROFESSIONAL TROLLEY



Pull-out tray for keyboard and mouse



Articulated arm for auxiliary monitor



Extendable tray with drawer

- 4 trays; 3 of them are height adjustable in steps of 25 mm.
- Strong shelves: 50 x 50 cm. Maximum load capacity: 50 kg
- Monitor mounting system VESA. 75 & 100.
- Pull-out tray for keyboard and mouse.
- 4 swivel wheels Ø 13 cm, two of them with brakes.
- Groove for cables.
- Dimensions: 162 x 66 x 50 cm.
- Weight: 75 kg.
- Colour: light grey.
- Options:
 - Independent isolation transformer.
 - Monitor articulated arm.
 - Extendable tray with drawer.
 - Holders for 3 rigid endoscopes, videocamera and fiber optic cable.
 - Support for flexible endoscope.



Specifications are subject to change without notice

